




PARASURE

ParaScope Protein Test Kit

For the rapid detection of protein residues in endoscope

- Conforms to recommendations in EN ISO 15883
- Results within 10 seconds
- Easy to use - no incubation required
- Clear colour change
- Detects protein residue inside endoscope channels
- 2.5 metre long swabs available with 3 different head sizes: 1.9mm, 2.8mm and 3.7mm



The ParaScope Protein Test Kit has been designed to allow the quick and reliable detection of residual proteins left behind inside endoscope channels. The long swabs provided are simply pushed down the length of the channel, removed and then swirled for 10 seconds in the liquid, to provide a clear and almost instantaneous result.

Based on a dye-binding solution used in clinical chemistry, this test is sensitive enough to detect protein residues to within 1 microgram. The protein test liquid will change to various shades of blue depending on the amount of protein detected, as can be seen on the colour guide overleaf. Any shade of blue indicates that protein has been detected.

For testing the surfaces of washer disinfectors or surgical instruments, use the ParaSure Protein Test, PS1700.

Order Details

PS1710

ParaScope 1.9mm Protein Test Kit
box of 25 tests with 25 2.5 metre long swabs
with 1.9mm head

PS1720

ParaScope 2.8mm Protein Test Kit
box of 25 tests with 25 2.5 metre long
swabs
with 2.8mm head

PS1720

ParaScope 3.7mm Protein Test Kit
box of 25 tests with 25 2.5 metre long
swabs
with 3.7mm head

Registered Address:

ParaSure Limited Units B10/11 Oakendene
Estate Bolney Rd Cowfold West Sussex RH13
8AZ England

+44 (0)1403 864777

Ref: ParaSure Tech Sheets 02/13 Rev.D



WWW.PARASURE.CO.UK



WWW.PARASUREORGANICS.COM



PARASURE

Multi Lumen Channel Control


TECHNICAL SPECIFICATIONS

Product name:	Protein Test Kit
Unit of measure:	box of 25 tests with 25 short swabs and 2 positive controls
Conforms to:	EN ISO 15883
Storage conditions:	if refrigerated, shelf life is extended from 6 months to 2 years
Frequency of use	as required
Number of uses:	single use
Supplied with:	instructions for use and record labels

Instructions for Use

1. If the object to be tested is already wet, vigorously swab the object, focusing in hinges or crevices which may be contaminated. If the object of interest is dry, simply wet the swab with tap water before vigorously swabbing.
2. Unscrew the cap and swirl the swab in the brown reagent for about 10 seconds, then discard the swab and screw the cap back on the reagent vial.
3. Visually inspect the reagent for colour change. If the reagent has turned blue, protein residue is present. The darker the blue colour the more protein has been detected. If the reagent remains brown, protein residue has not been detected.

Two test strips have been supplied for use as a positive control. To use them, simply wet the swab, and vigorously rub it on the red test soil on the aluminium strip. Unscrew the cap and swirl the swab in the reagent. Discard the swab and screw the cap back on the reagent vial. The result will be bright blue, and provides a reference colour for a highly positive reaction to protein.



For further information about any of our products, please visit:
www.parasure.co.uk



WWW.PARASURE.CO.UK



WWW.PARASUREORGANICS.COM

Registered Address:

ParaSure Limited Units B10/11 Oakendene
Estate Bolney Rd Cowfold West Sussex RH13
8AZ England

+44 (0)1403 864777