



# PARASURE

## Seal Sure

- Routine monitor for sealing device
- In accordance with ENISO 11607-2 2006
- Records sealing time, temperature and time
- Cost effective
- Easy to interpret



## Seal Shield

The ParaSure Seal Sure is for use as a routine monitor of a sealing device and for operational qualification within the framework of process validation in accordance with ENISO 11607-2 2006.

It is used to ensure optimal sealing seams are achieved by the heat sealing machine by monitoring sealing time, temperature and pressure.

In accordance with EN ISO 11607-2 the sealing seams must show the following.

- Intact sealing for a specified seal width
- No channels or open seals
- No punctures or tears
- No material delamination or separation.

The ParaSure Seal Sure ensures that the above properties are fulfilled and therefore a subjective (i.e. purely visible) interpretation of the sealing seam is avoided and human error excluded.

The test should be performed daily at the beginning of the working day and recorded as part of the CSSD record systems for the entire decontamination process.

### Order Details

#### PS1900

ParaSure  
Face Shields

Box of 200



[WWW.PARASURE.CO.UK](http://WWW.PARASURE.CO.UK)



[WWW.PARASUREORGANICS.COM](http://WWW.PARASUREORGANICS.COM)

#### Registered Address:

ParaSure Limited Units B10/11 Oakendene  
Estate Bolney Rd Cowfold West Sussex RH13  
8AZ England

+44 (0)1403 864777

Ref: ParaSure Tech Sheets 02/13 Rev.D



# PARASURE

## Seal Sure

### TECHNICAL SPECIFICATIONS

<b>Name:</b>	Seal Sure
<b>Conformity:</b>	EN-ISO 11607-2 2006
<b>Correct Symbols</b>	In accordance with EI-ISO 11607-2 2006
<b>Disp/single use:</b>	yes
<b>Extra Information:</b>	
<b>Non-toxic/Lead free</b>	yes
<b>Recordable Information:</b>	
<b>Sealing Machine</b>	yes
<b>Location</b>	yes
<b>For approval</b>	yes
<b>Date</b>	yes
<b>Remarks</b>	yes
<b>Shelf Life</b>	5 years



For further information about any of our products, please visit:  
[www.parasure.co.uk](http://www.parasure.co.uk)



[WWW.PARASURE.CO.UK](http://WWW.PARASURE.CO.UK)



[WWW.PARASUREORGANICS.COM](http://WWW.PARASUREORGANICS.COM)

### Registered Address:

ParaSure Limited Units B10/11 Oakendene  
Estate Bolney Rd Cowfold West Sussex RH13  
8AZ England

+44 (0)1403 864777

Ref: Parasure Tech Sheets 02/13 Rev.D