



# PARASURE

## Type 6 Self-Adhesive Indicator

18 minutes at 134°C

- Conforms to ISO 11140-1:2014
- Specifically calibrated to 3.5 minutes at 134°C
- Monitors Time, Steam & Temperature
- Easy to interpret results
- Small, compact design
- Self-adhesive - easy to record



The ParaSure Type 6 emulating chemical indicators offer the highest level of proof of the achievement of sterility assurance and are specifically calibrated to the cycle requirements of 134°C for 3.5 minutes.

The indicators offer validated security that every tray or item has been effectively processed, thereby providing a direct link between the sterility of the items and the patient.

A clear colour change and self-adhesive backing ensure results which are easy to interpret and to record.

### Order Details

#### PS1210

Type 6 Self-Adhesive Indicators

bag of 250 indicators

### Registered Address:

ParaSure Limited Units B10/11 Oakendene  
Estate Bolney Rd Cowfold West Sussex RH13  
8AZ England

+44 (0)1403 864777



[WWW.PARASURE.CO.UK](http://WWW.PARASURE.CO.UK)



[WWW.PARASUREORGANICS.COM](http://WWW.PARASUREORGANICS.COM)



# PARASURE

---

## Type 6 Self-Adhesive Indicator

3.5 minutes at 134°C

### Technical Specifications

<b>Product name:</b>	Type 6 Self-Adhesive Indicator
<b>Unit of measure:</b>	1 bag of 250 indicators
<b>Conforms to:</b>	EN ISO 11140-1
<b>Frequency of use:</b>	daily
<b>Number of uses:</b>	single use
<b>LOT number:</b>	on bag and indicator
<b>Expiry date:</b>	on bag and indicator
<b>Supplied with:</b>	instructions for Use

For further information about any of our products, please visit:  
[www.parasure.co.uk](http://www.parasure.co.uk)

---



[WWW.PARASURE.CO.UK](http://WWW.PARASURE.CO.UK)



[WWW.PARASUREORGANICS.COM](http://WWW.PARASUREORGANICS.COM)

#### Registered Address:

ParaSure Limited Units B10/11 Oakendene  
Estate Bolney Rd Cowfold West Sussex RH13  
8AZ England

+44 (0)1403 864777