

## **Autoclave Tape**

- Conforms to ISO 11140-1 Class 1
- Easy to interpret
- Cost effective
- Individual rolls marked with lot number and expiry date
- Non-toxic
- Lead free



### **Autoclave Tape**

- Autoclave tape for steam sterilisation enables the user to see at a glance that an item has been processed. The easy to interpret black stripes prove that the pack has been exposed to the steam sterilisation process.

The tape securely sticks to non woven, paper, board, metal, plastic and glass and is manufactured from treated crepe paper and coated with high performance cross linked heat sensitive adhesive

# Order Details

#### PS1510

19 mm Autoclave tape Cartons of 48 rolls

#### PS1520

25 mm Autoclave tape Cartons of 36 rolls





#### Name: Autoclave tape for steam Sterilisation Classification: EN-ISO 11140-1 (Class 1) Operating temp: 121°C - 134°C As per specific sterilisation cycles minimum Time: 134°C 3.5 minutes or equivalent Storage conditions: 10°C - 30°C RH: 30 - 70% Lot no. yes, on the core Prod date: yes, on the outer carton Exp date: yes, on the core Correct symbols: in accordance with the standard Indicator: yes **Expiry Date:** 2 years **Extra information:** Medical crepe paper Substrate Coated Information can be written on the tape with a Recordable Information: ParaSure PIP Colour change: white to black 25mm width tape = bags of 36 rolls Packaging: 19mm width tape = bags of 48 rolls

## **Autoclave Tape**



Address: Squirrels | Kithurst Park | Storrington | West Sussex | RH20 4JH | England | Tel: +44 (0) 1903 741390



