




## STOP! be sure...be ParaSure

### Type 6 Self-Adhesive Indicators



**Unprocessed indicator**  
Ready to use



**PASS**

Items in tray/pack can be safely used

**PASS**

Items in tray/pack can be safely used

**Successfully processed indicator**  
Colour has changed to teal or darker



**FAIL**  
Do NOT use items in tray/pack

**Unsuccessfully processed indicator**  
Evidence of original pink colour



**FAIL**  
Do NOT use items in tray/pack

**Unsuccessfully processed indicator**  
Evidence of original pink colour



**FAIL**  
Do NOT use items in tray/pack

**Unsuccessfully processed indicator**  
Evidence of original pink colour

PS1210 3.5 minutes at 134 degrees	PS1230 5.3 minutes at 134 degrees
PS1220 4.0 minutes at 134 degrees	PS1240 7.0 minutes at 134 degrees

