

Ultrasonic Foil Test

To be used with PS 2022



Ultrasonic foil tests consist of foil-based material. If the appliance operates correctly, the cavitation energy emitted from sources will dismember the foil inside the frame.

The wash test should be used daily and the results should be recorded. There are areas on the product that can be filled in so that the results can be recorded.

PS 2220 50 pieces per box

Order Details

W W W . P A R A S U R E . C O . U K

V W W . P A R A S U R E O R G A N I C S . C O M

Registered Address:

ParaSure Limited Units B10/11 Oakendene Estate Bolney Rd Cowfold West Sussex RH13 8AZ England

> +44 (0)1403 864777 Ref: ParaSure Tech Sheets 12/12 Rev.C




PARASURE

Ultrasonic Foil Test

To be used in conjunction with PS 2022

TECHNICAL SPECIFICATIONS

Physical Parameters	
Process	Ultrasonic Washer
Packaging	50 PCS. / BOX
Speciality	Initial Situation: Full Foil Signal Situation: Fragmented Foil



For further information about any of our products, please visit: www.parasure.co.uk



W W W . P A R A S U R E . C O . U K

Registered Address:

ParaSure Limited Units B10/11 Oakendene Estate Bolney Rd Cowfold West Sussex RH13 8AZ England

+44 (0)1403 864777