




PARASURE

Plasma Sterilisation Monitoring Indicator

- Calibrated for VHP sterilisation processes
- Conforms to EN ISO 11140-1 Type 4
- Non-toxic
- Clear colour change
- Easy to interpret - easy to record
- Lead free



Plasma Sterilisation Monitoring Indicator

The ParaSure Plasma Sterilisation Monitoring indicator is designed for use in a VHP steriliser.

The indicator will exhibit a clear colour change when exposed to hydrogen peroxide during the diffusion stage of the process. When utilised properly, the VHP indicator will give clear visual confirmation that the item was exposed to hydrogen peroxide

Order Details

PS1800

Plasma
Sterilisation
Monitoring
Indicator
Box of 250



WWW.PARASURE.CO.UK



WWW.PARASUREORGANICS.COM

Registered Address:

ParaSure Limited Units B10/11 Oakendene
Estate Bolney Rd Cowfold West Sussex RH13
8AZ England

+44 (0)1403 864777

Ref: ParaSure Tech Sheets 08/12 Rev.B



PARASURE

Plasma Sterilisation Monitoring Indicator

TECHNICAL SPECIFICATIONS

Name:	Plasma Sterilisation Monitoring Indicator
Classification:	N/A FDA cleared
Time:	30 minutes
Storage conditions:	50 - 100°F (10 - 38°C)
RH:	10% - 90%
Lot no.	yes
Prod date:	yes on box
Exp date:	2 years
Indicator:	yes, changes from blue to pink
Disp/single use:	yes
Extra Information: Non-toxic/Lead free Indicator Substrate material	yes polystyrene
Recordable Information: Date: Operator Pass/fail result Colour change	yes yes yes blue to pink



Before
Original Indicator



After
After total process

For further information about any of our products, please visit:
www.parasure.co.uk



WWW.PARASURE.CO.UK



WWW.PARASUREORGANICS.COM

Registered Address:

ParaSure Limited Units B10/11 Oakendene
Estate Bolney Rd Cowfold West Sussex RH13
8AZ England

+44 (0)1403 864777